

Figure 1




Figure 2

FIG. 3




Figure 3

treatment
vessels




Fig. 4A




Fig. 4B




Fig. 4C

pre-treatment
assembly




Fig. 5A




Fig. 5B




Fig. 5C

post-treatment
transfer




Fig. 6A




Fig. 6B

vacuum extract fluid and tissue




Figure 7




Figure 8




Figure 9

Figure 11

SYSTEM SPECIFICATIONS	EXTRACTION	TRANSFORMATION	RESEARCH
INSTRUMENT CONTROL: LabVIEW			
x-y-z positioning (sample)	yes	yes	yes
z' axis (transducer)	manual, 25mm range	manual, optional auto	manual
Temperature feedback to protocol	yes	yes	yes
partial treatments	yes	optional	no
cavitation detection	no	optional	yes
video detection and analysis	no	optional	yes
USER INTERFACE: LabVIEW			
treatment protocol	fixed	user adjustable	flexible
select treatment positions	pre addressed	user adjustable	flexible
temperature profile record	optional	yes	yes
timing information	yes	yes	yes
ELECTRICAL:			
Power: 110V, 20A			
EQUIPMENT:			
Chiller	yes	no	yes
RF Amplifier	yes	yes	yes
Arbitrary waveform generator	yes	yes	yes
oscilloscope	no	optional	yes
Computer	yes	yes	yes
motion control	yes	yes	yes
I/o boards			
amplifier	yes	yes	yes
xy stage	no	optional	yes
IR temperature measurement	yes	yes	yes
video	no	optional	yes
laser sight/cross-hairs	yes	yes	yes
vacuum fixture	yes	no	no
Transducer			
matching network			
cables			
circulation pump			
convection cooling			
filter			
cavitation detection	no	optional	yes

Figure 12

LabVIEW PROGRAMMING TASKS

GENERAL

	Extraction	Transformation
display revision level	x	x
safety interlocks	x	x
time and date stamp		x
STOP function	x	x
save configuration to file	user can reset defaults	x
operating parameters		x ^f
protocol		x
save data to file		
treatment positions and protocols		x
temperature profile		x
error conditions		x
password protection on Vis	x	x
load configuration from file		x
user selects treatment positions	x	x

DISPLAY

User selectable treatment positions -graphical	x	x
current status		
treatment position -graphical	x	x
current protocol	by name	x
-voltage		x
-duty cycle		x
-etc		x
time to finish current sample	x	x
safety interlock status	x	x
sample temperature, graph and current temp		x
time and date		x

ULTRASONICS

initialize instrument(s)	x	x
stop function	x	x
mix and treat	predetermined	userprogrammable
frequency	predetermined	x
voltage-treat	predetermined	x
voltage-mix		x
pulselength-treat	predetermined	x
pulselength-mix		x
deadtime-mix>treat		x
deadtime-treat>mix		x
Total cycles (or time)	predetermined	x
cavitation detection		optional

POSITIONING

setup and diagnostics		
initialize stepper control board	x	x
calibrate (home)	x	x
check limits (limit switches)	x	x

PCT/US99/25274

Transformation

predetermined
x

x
x
x
x

go to next well at set temperature rise

Figure 13